



# SimpliMet™ 4000 Service Plans



## Protect Your Investment

Labs run more efficiently and consistently with well-maintained equipment.



## Reduce Unplanned Downtime

Regular equipment maintenance will reduce the hassle of unexpected downtime and associated costs.



## Ensure Consistent Results

Confidently analyze lab results with the equipment maintained by the trusted experts at Buehler.



## Service Plan Options

Identify the best Service Plan based on the lab's needs.

	Preventative Maintenance	Extended Warranty
Visits Per Year	1	2
Wear Parts Kit	✓	✓
Unlimited Phone Support	✓	✓
Replacement Part Discount	10%	25%
Unlimited Labor Repairs		✓
Priority Response		✓

- Preventative Maintenance ensures consistent results
- Visit includes Wear Parts Kit and 14 Point Maintenance Check list *(see reverse for more details)*
- Extended Warranty must be purchased prior to the expiration of active warranty



# Maintenance Checklist

Our comprehensive 14 Point Maintenance Checklist is thorough, ensuring the SimpliMet 4000 is running efficiently.

## 14 Point Maintenance Checklist

- ☑ Check and inspect water & drain hoses
- ☑ Check water solenoid for proper operation
- ☑ Check inside the unit for oil leaks
- ☑ Vacuum inside the unit and clean all surfaces
- ☑ Check bayonet closure assembly for proper operation
- ☑ Check inside the unit for any water leaks
- ☑ Check PCB boards, ribbon cables and wire terminals
- ☑ Clean top & bottom rams.
- ☑ Check cooling fan for proper operation
- ☑ Inspect molds to verify unit is working properly
- ☑ Record hours, temperature and voltage readings
- ☑ Check for error codes
- ☑ Check and update firmware version
- ☑ Check all functions of machine for proper operation



## Wear Parts Kit

The parts listed below are included with every Preventative Maintenance visit

**SimpliMet 4000 Wear Parts Kit**      201011SRV

### *Includes*

- Water Solenoid
- Water In Quick Disconnect Plug
- Water Out Quick Disconnect Plug
- Hose Clamp
- Mold Bolt



**BUEHLER Worldwide Headquarters**  
**North America-South America Offices**  
 41 Waukegan Road  
 Lake Bluff, Illinois 60044-1699 USA  
 P: 800 BUEHLER (800-283-4537)  
 P: (847) 295-6500  
 www.buehler.com | info@buehler.com

**European Headquarters**  
**BUEHLER Germany**  
 info.eu@buehler.com

**BUEHLER France**  
 info.eu@buehler.com

**BUEHLER United Kingdom**  
 info.eu@buehler.com

**BUEHLER China**  
 info.cn@buehler.com

**BUEHLER Japan**  
 info.japan@buehler.com

**BUEHLER Asia-Pacific**  
 info.asia@buehler.com